



BUCKY TABLE

BK-12HK



Features**1**

The table top slides longitudinally and laterally in a wide range. Positioning can be easily carried out. Also the table slides only longitudinally when using a button switch.

2

Cassettes from 14"×17" to 6¹/₂"×8¹/₂" in size are usable for both horizontal and vertical positions.


*) Cassettes are available in cms or inches in accordance with the IEC and also have a thickness of 15 mm or less

OPTIONS

- Oppression band
- CFRP tabletop (810×2200mm)
- Side cassette holder
- Handle
- Photo timer (Automatic exposure control devices)

SPECIFICATIONS

Size of table mm (inch)	810×2200 (31.9×86.6)
Material of tabletop	WOOD
Table height from floor mm (inch)	700 (27.6)
Tabletop travel mm (inch)	
Longitudinal	1100 (43.3)
Lateral	±125 (±4.9)
Tabletop lock	Electromagnetic lock (It is released when activated)
Bucky stroke mm (inch)	350 (13.8)
Size of loadable cassettes	14"×17" - 6 ¹ / ₂ "-8 ¹ / ₂ " usable for both horizontal and vertical positions
Bucky unit	Grid ratio R=12:1
	Grid density N=40lines/cm
	Focal distance fo=100cm
	Spacer AL
	Size 479×438cm
Power requirements (supplied from an x-ray generator)	Single phase 100V, 1A, 50/60Hz
Weight	120kg
Maximum allowable load	200kg

DIMENSIONS

unit:mm

**SHIMADZU CORPORATION. International Marketing Division**

3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan Phone: 81(3)3219-5641 Fax: 81(3)3219-5710

SHIMADZU MEDICAL SYSTEMS

20101 South Vermont Ave., Torrance, California 90502, U.S.A.
Phone: 1(310)217-8855 Fax: 1(310)217-0661

SHIMADZU DEUTSCHLAND GmbH

Albert-Hahn-Strasse 6-10, D-47269 Duisburg, F.R. Germany Phone: 49(203)7687-0 Fax: 49(203)766625

SHIMADZU (ASIA PACIFIC) PTE LTD.

16 Science Park Drive #01-01 Singapore Science Park, Singapore 118227, Republic of Singapore
Phone: 65-6778 6280 Fax: 65-6779 2935

SHIMADZU (HONG KONG) LTD.

Suite 1028 Ocean Center, Harbour City, Tsim Sha Tsui, Kowloon, Hong Kong
Phone: 852-2375-4979 Fax: 852-2199-7438

SHIMADZU INTERNATIONAL TRADING (SHANGHAI) Co., LTD. SHANGHAI OFFICE

24th Floor, Shanghai Xin Hualian Building, No.755 Huaihai Zhong Lu, Shanghai, China
Phone: 86-21-6472-8442 Fax: 86-21-6472-8648

SHIMADZU MEDICAL SYSTEMS (OCEANIA) PTY. LTD.

Units E, 10-16 South Street, Rydalmere N.S.W. 2116, Australia
Phone: 61(2)9898-2444 Fax: 61(2)9684-4055


SHIMADZU DO BRASIL COMÉRCIO LTDA.

Avenida Marquês de São Vicente, 1771. Barra Funda CEP : 01139-003-São Paulo-SP, Brasil
Phone: (55)11-3611-1688 Fax: (55)11-3611-2209

Overseas Offices

Istanbul, Moscow

URL <http://www.shimadzu.com>



Shimadzu Corporation Medical Systems Group has been certified by TÜV Rheinland as a manufacturer of medical equipment and systems in compliance with ISO9001:2000 Quality Management Systems and ISO13485:2003 Medical Equipment Quality Management Systems.

Remarks;

- *Every value in this catalogue is a standard value, and it may vary a little from the actual at each site.
- *The appearances and specifications are subject to change for reasons of improvement without notice.
- *Certain configurations may not be available pending regulatory clearance. Contact your Shimadzu representative for information on specific configurations.
- *Before operating this system, you should first thoroughly review the Instruction Manual.